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
## Polymorphic Transitions in Carbamazepine by Low Frequency probe-type Raman spectroscopy

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
Fig. 1

Table 1

© 2019 COHERENT COMPANY. All Rights Reserved.

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
Fig. 6

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- (1) Láng, P., Kiss, V., Ambrus, R., Farkas, G., Szabó-Révész, P., Aigner, Z., & Várkonyi, E. (2013). Polymorph screening of an active material. *Journal of pharmaceutical and biomedical analysis*, 84, 177-183.
- (2) 今井昭生. (2019). 合成原薬の医薬品ライフサイクルマネジメントの実際—アルツハイマー型認知症治療薬アリセプト原薬（ドネペジル塩酸塩）の事例—. *日本PDA学術誌GMPとバリデーション*, 21(1), 1-21.
- (3) 小針昌則 Masanori 日揮技術ジャーナル Vol.3 No.4 (2014) 晶析プロセスの新たな解釈と溶媒媒介転移
- (4) Starbuck, C., Spartalis, A., Wai, L., Wang, J., Fernandez, P., Lindemann, C. M., ... & Ge, Z. (2002). Process optimization of a complex pharmaceutical polymorphic system via in situ Raman spectroscopy. *Crystal Growth & Design*, 2(6), 515-522.
- (5) Inoue, M., Hisada, H., Takatori, K., Koide, T., Fukami, T., Roy, A., & Carriere, J. (2020). Solid-state analysis of alpha-cyclodextrin inclusion complexes using low-frequency Raman spectroscopy. *Analytical Chemistry*, 93(2), 704-708.
- (6) Hédoux, Alain, et al. Polymorphic transformation of anhydrous caffeine upon grinding and hydrostatic pressurizing analyzed by low-frequency raman spectroscopy, *Journal of pharmaceutical sciences*, 102, 1, 162-170 (2013).
- (7) Hisada, Hiroshi, et al. Determining the Distribution of Active Pharmaceutical Ingredients in Combination Tablets Using Near IR and Low-Frequency Raman Spectroscopy Imaging, *Chemical and Pharmaceutical Bulletin*, 68, 2, 155-160 (2020).
- (8) Remoto, P. I. J., Bérziš, K., Fraser-Miller, S. J., Korter, T. M., Rades, T., Rantanen, J., & Gordon, K. C. (2023). Exploring the Solid-State Landscape of Carbamazepine during Dehydration: A Low Frequency Raman Spectroscopy Perspective. *Pharmaceutics*, 15(5), 1526.
- (9) O'Brien, L. E., Timmins, P., Williams, A. C., & York, P. (2004). Use of in situ FT-Raman spectroscopy to study the kinetics of the transformation of carbamazepine polymorphs. *Journal of pharmaceutical and biomedical analysis*, 36(2), 335-340.
- (10) Hisada, H., Okayama, A., Hoshino, T., Carriere, J., Koide, T., Yamamoto, Y., & Fukami, T. (2020). Determining the distribution of active pharmaceutical ingredients in combination tablets using near IR and low-frequency Raman spectroscopy imaging. *Chemical and Pharmaceutical Bulletin*, 68(2), 155-160.

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© 2019 MarqMetrix. All Rights Reserved.